

## CRK-7000 AUTO REF-KERATOMETER CR-7000 AUTO REFRACTOMETER



CRK-7000



CR-7000

**Huvitz**  
Pacing Progress toward People

US VISION ON TECHNOLOGY  
**Jaggi**

The Charops CR / CRK 7000 is the latest innovation from the parent organization of M/s Huvitz Co Ltd, Korea. The Brand Charops is being re-launched by M/s Huvitz Co Ltd, Korea to promote economical, feature-rich and sturdy machines. The Equipment have been manufactured under the strict quality control of M/s Huvitz Co Ltd, Korea to give you the most feature ready range of instruments at the most affordable prices.

CRK-7000



**WIDER MEASURING RANGE**  
Superior performance with a wider measuring range in Refraction & Keratometry readings is provided.



**PUPIL & IRIS SIZE MEASUREMENT**  
Pupillary and Corneal Size can be measured for precise fitting of Contact Lens and help in refractive surgery.


CR-7000



**5.7 INCH COLOUR LCD**  
Displays sharper & clearer image on screen on a real time basis.



**CONVENIENT ONE TOUCH LOCK**  
Locking & unlocking is done easily through one touch



**FOCUSING INDICATOR**  
Focus for accurate reading shows up when properly centered, thus providing accurate results



**HIGH SPEED PRINTER**  
Thermal High speed printer provides with fast print-outs.

### SPECIFICATIONS

	CR 7000 Autorefractometer	CRK 7000 Autorefractometer
Vertex Distance	0,12,13,5,15 mm	
Spherical value	-25.00 Ds to +22.00 Ds	
Cylinder value / axis	Upto +/- 10.00 Dc, axis upto 180 in 1deg steps	
Pupillary Distance	10 - 85 mm in 1mm steps	
Minimum Pupil Diameter	2,00 mm	
Ref Mode	Available	
Ker Mode	Not available	Available
Ref / Ke r mode	Not available	Available
CLBC mode	Not available	Available
Ker range	NA	33-67.50 D
Display mode of data	Available for the last 10 readings taken of each eye	
Internal Printer	Available	
Size mode for HVID	Available	
Display	5.7 inch Color LCD monitor	
Power Supply	AC 100 - 220v ac , 50/60 Hz	
Dimensions	275(W) X 510(D) X 450 (H) / 20kgs	
Power Saving Mode	Available , auto off in 5 minutes	

Manufactured by  
**Huvitz**  
Pacing Progress toward People  
Huvitz Building 488-3, Geumjeong-dong  
Guropo-si, Gyeonggi-do, 435-862, Korea  
Tel : 82-31-442-8888 Fax : 82-31-477-8617  
http : //www.huvitz.com

Exclusively Marketed in India by :

US VISION ON TECHNOLOGY  
**Jaggi** JAGGI BROTHERS  
311, Pratap Chambers, Gurukulam Road, Karol Bagh, New Delhi - 110005  
Tel : +91-11-45052669, 28757334, 45052612  
E-mail : sales@jaggi.com Website : www.jaggi.com