

Milan Optical Devices



AR - 800 A / ARK - 800
Auto Ref / Keratometer





AR - 800 A / ARK - 800 Auto Ref / Keratometer



FEATURES

- ❖ IT ADOPTS EMBEDDED MICRO-PROCESS AND FAST IMAGE PROCESS TECHNIQUE, WHICH TAKE MEASUREMENT QUICKLY
- ❖ WITH AUTO – READING FEATURES , OPERATOR CAN TAKE MEASUREMENT AUTOMATICALLY ONCE EQUIPMENT OBTAINS SATISFACTORY FOCUS AND ALIGNMENT
- ❖ QUICK MEASUREMENT OF DIOPTER FOR MINIMUM PUPIL SIZE
- ❖ NEW DESIGN OF COLOURFULL SCREEN (5.7' TFT LCD)
- ❖ AUTO FOGGING CHART
- ❖ PUPIL DISTANCE AUTO MEASUREMENT AVAILABLE
- ❖ CHOICE FOR VERTEX DISTANCE

MEASURING RANGE	SPHERICAL	-20 ~ +20D (VD = 12) 0.12 STEPS
	CYLINDER	-8 ~ +8D 0.12D STEPS
	AXIS	0 ~ 180

PUPIL DISTANCE	45 ~ 85 MM , 1 MM STEPS
-----------------------	--------------------------------

VERTEX DISTANCE	0 MM , 12 MM , 13.75 MM
------------------------	--------------------------------

MIN PUPIL SIZE	2.2 MM
-----------------------	---------------

CHART	AUTO FOG
--------------	-----------------

DISPLAY	5.7 ' TFT LCD
----------------	----------------------

PRINT	THERMAL PRINTER
--------------	------------------------

POWER SUPPLY	100 ~ 120 V / 200 – 240 V , 50/60 HZ
---------------------	---

DIMENSION / WEIGHT	480 (L) X 282 (W) X 445 (H) MM / 17 KG
---------------------------	---